


Scholten

Bioptome

Phone: Toll Free 866 452 3525 • Fax: (209) 365 1437 • www.bioptome.com • 170 Commerce Street , Suite 101, Lodi , California 95240 USA



Reusable Endomyocardial Biopsy Forceps

<u>Catalog Number</u>	<u>French Size</u>	<u>Length</u>	<u>Curved Tip</u>	<u>Jaw Cup Dia.</u>
<u>101</u>	<u>9</u>	<u>50cm</u>	<u>Yes</u>	<u>2.8mm</u>
<u>102</u>	<u>8</u>	<u>50cm</u>	<u>Yes</u>	<u>2.5mm</u>
<u>103</u>	<u>6.5</u>	<u>50cm</u>	<u>Yes</u>	<u>2.1mm</u>
<u>105</u>	<u>6.5</u>	<u>100cm</u>	<u>No</u>	<u>2.1mm</u>
<u>115</u>	<u>9</u>	<u>100cm</u>	<u>No</u>	<u>2.8mm</u>
<u>116</u>	<u>8</u>	<u>100cm</u>	<u>No</u>	<u>2.5mm</u>
<u>106</u>	*	*	*	*


* Custom catheter lengths are available for either right or left ventricular biopsy forceps upon request at no extra cost

The reusable Scholten Bioptome is the high performance and cost effective approach to endomyocardial biopsy. As expressed by cardiologists and other health professionals, the Scholten Bioptome yields more and better quality specimens at a lower cost.

Features & Performance

- A design proven by over 30 years of successful clinical use.
- Industry standard cutting head diameters and shaft lengths are available.
- Manufactured using biocompatible materials.
- Unique hollow cut cutting edge for superior sharpness.
- Single action cutting jaw (single moving jaw) provides superior cutting force.
- Super precision machining provides precise cutting edge alignment for every cut.
- Familiar locking finger ring actuation handle.
- Formable distal curve.
- Optional Bioptome reconditioning service available through SSI.

From WA Baumgartner, Heart and Lung Transplantation, 2nd Edition, P.273; Copyright Elsevier, 2002.



Bioptome Factory Reconditioning Service

Precision manufacturing enables the reusable Scholten Bioptome to be reconditioned numerous times, resulting in many years of service. Reconditioning includes flexible shaft replacement, inspection and careful evaluation of each component. At a remarkably low cost, the instrument is guaranteed to perform like new when reconditioned.



Scholten Surgical Instruments, Inc.

ISO 13485 Certified Quality Management System

www.bioptome.com