

Bioptome Accessories

Phone: Toll Free 866 452 3525 • Fax: (209) 365 1437 • www.bioptome.com • 170 Commerce Street, Suite 101, Lodi, California 95240 USA

Bioptome Cleaning Device #2

The Bioptome Cleaning Device offers the advantages of ease of use and portability, and uses a standard syringe.



Product Codes:

108	Cleaning Device # 2 -Bioptome flex shaft cleaning device
-----	--

- Cleans and lubricates each catheter in a relatively short period of time.
- Simple to use. Comes complete with instructions.
- Prompt use enables Bioptomes to obtain the maximum number of biopsies before undergoing reconditioning.

Bioptome Autoclaving/Storage Plate

- Important for the protection of your instruments during sterilization and storage



Product Codes:

112	Autoclaving / Storage Plate for 50cm Length bioptome
113	Autoclaving / Storage Plate for 100 – 103cm Length Bioptome
114	Autoclaving / Storage Plate for custom size bioptome



Scholten Surgical Instruments, Inc.

ISO 13485 Certified Quality Management System

www.bioptome.com