


Scholten

Novatome™

Single-use endomyocardial biopsy forceps

A biopsy design proven by over 30 years of successful use. The same sharp and steerable advantages in a new disposable format.




Features and Performance

- Perfect balance of flexibility and stiffness offering precise navigation and positioning.
- Formable distal curve.
- Familiar locking finger ring actuation handle.
- Super precision machining provides precise cutting edge alignment for every cut.
- Single action cutting jaw (Single moving jaw) provides superior cutting force.
- Manufactured using biocompatible materials.
- Sterile and disposable for maximum patient safety and facility convenience without compromising performance and strength.
- Packaged sterile in protective tray with easy peel pouch.
- Industry standard cutting head diameters and shaft lengths available.

Catalog Number	French size	Length	Curved tip	Jaw cup dia.
SU101-50	9	50cm	X	2.84 mm
SU102-50	8	50cm	X	2.51 mm
SU103-50 *	7	50cm	X	2.28 mm
SU104-50 *	6	50cm	X	1.99 mm
SU101-100 *	9	100cm		2.84 mm
SU102-100 *	8	100cm		2.51 mm
SU103-100 *	7	100cm		2.28 mm
SU104-100 *	6	100cm		1.99 mm

* Available Soon



Tradition meets the Future

In the late seventies, the Stanford Heart Team asked Jake Scholten to provide them with a high quality endomyocardial biopsy forceps. That instrument, known as the **Scholten Bioptome**, has been trusted for its uncompromising quality and reliability. The latest Scholten disposable Bioptome model is called the **Novatome™** and, like its predecessors, is specially designed to acquire maximum specimen size for more accurate pathology results. With a reputation for its ease of positioning and heart rhythm feedback, Scholten Bioptomes have earned high praise.

The new Novatome™ disposable model makes Scholten the best choice for safety and precision. Clearly, with a pedigree of more than 30 years of clinical experience, the new disposable model Novatome™ will continue to set the standard worldwide.

Scholten Surgical Instruments, Inc. • ISO 13485 Certified Quality Management

Phone: Toll Free 866 452 3525 • Fax: (209) 365 1437 • www.bioptome.com • 170 Commerce Street, Suite 101, Lodi, California 95240 USA